



GS-Series Leather Drum Dyeing (lab use)



Note :
Due to keeping innovation,
the specification and
Machine pictures or color
are subject to alteration
without notice!

Application :

GS-series Leather Drum Dyeing (lab use) is an essential experiment equipment for modern leather production, and mainly used for small batch and multi-pieces test of various kinds of leather to obtain the best production processes, such as preparation, tanning, dyeing etc.

Parameters :

Model	Unit	GS-D800	GS-D1000	GS-D1200	GS-D1400	GS-D1600	GS-D1800
Drum diameter	mm	800	1000	1200	1400	1600	1800
Length of drum	mm	400	500	600	600	600	600
Volume of drum	L	200	390	680	920	1200	1520
Loading capacity (leather cut)	KG	16	32	54	74	96	120
Drum speed	r/min	0-25	0-25	0-20	0-16	0-16	0-16
Motor power	KW	1.5	1.5	1.5	3	4	5.5
Heating power	KW	4.5	6	6	6	6	9
Temperature range	°C	10~100±2					
Power		Motor 50HZ 380V, Heater 50HZ 380V					
Dimension	cm	130x100x150	140x130x150	160x140x170	180x140x180	200x150x210	220x160x220
Total weight	T	0.30	0.60	0.70	0.80	1.20	1.50