

© 2004 Blackwell Publishing Ltd *Journal of Internal Medicine* 255: 105–112



0000-0000-0000-0000

Items	Units	8/8/2020	4/8/2020
Max product volume	L	300	1000 + 3000
Discharge	PCU	500	500
Power demand	kW	100	100+0
Power L/D ratio	L/D	50	50
Disch. (Horsepower)	HP	100	100+0
Disch. heating capacity	HP	50	50+0
No. of heating coils	COIL	8	8+0
Coils of 2000	COIL	500	500+500
CR pump drive power	HP	20	20
Discharge flow	g/s	500	500

Optimal function for specific time periods such as: Morning/day, Day/evening, Weekend support, time and date some travel advice, will remain open until then, additional advice, new travel information



Due to keeping innovation
the specification and
Machine pictures or color
are subject to alteration
without notice!

Param.	Unit	WSD 20	WSD 115M
Max. dia of hole	mm	1200-1400	1800-2000
Maximal string	L	20-30	40-100
Template size	mm	1000-1400	1800-2000
Max. stroke	mm	800-1100	800-1000
Max. dia. diameter	mm	700	1000
Che. feeding power	kw	30	60
Che. feeding line segments	mm	8	8
String pressure	gpa	0.8	1.2
Air volume	m ³ /min	1.8	8
Cooling water pressure	MPa	0.3	0.3
Water consumption	L/min	180	240
Size of machine	LxWxH	2300x1500x3	1100x750x3
Weight of machine	TON	30	75